

Logi™ Range 3mm

Resposable™ Laparoscopic Instrumentation

Logi™ Range 3mm Instruments are a Resposable™ high quality, cost-effective solution for ultra-minimally invasive surgery.

The aim of 3mm surgery is to substantially reduce post-operative pain and trauma for the patient, promote faster recovery time and improve 'cosmesis' (minimal scarring) by making smaller incisions.

Key Strengths of Logi™ Range 3mm

- + Resposable™ - reusable handles and single use inserts
- + Fits our Logic™ reusable handles, including the 'Pistol-Grip' handle and patented 'Horizontal' handle both available with or without ratchet
- + Impressive strength despite reduced size
- + Perfect complement to YelloPort+plus™ 3mm laparoscopic port access system

Logi™ Grasp

Logi™ Dissect

Logi™
Extended Grasp

Logic™ Pistol Grip Handle
with Ratchet

For further detail about Logi™ Range 3mm please contact:

sales@surginno.co.uk

+44(0)113 230 7597

YelloPort+plus™ 3mm

Resposable™ Laparoscopic 3mm Port Access System

YelloPort+plus™ 3mm is a Resposable™ high quality, cost-effective laparoscopic port access system designed for use in ultra-minimally invasive surgery.

The aim of 3mm surgery is to substantially reduce post-operative pain and trauma for the patient, promote faster recovery time and improve 'cosmesis' (minimal scarring) through smaller incisions.



Key Strengths of YelloPort+plus™ 3mm

- + Resposable™ - reusable Cannulas and Obturators, and disposable sterile, pre-lubricated valves
- + Impressive strength despite reduced size
- + Selection of Cannulas and Obturators available
- + Perfect complement to 3mm Logi™ Range Instruments

For further detail about YelloPort+plus™ 3mm please contact:

sales@surginno.co.uk

+44(0)113 230 7597

PretzelFlex™

part of the EndoFlex™ range


Surgical Retraction System

PretzelFlex™ is the world's first reusable pretzel shaped organ and tissue retractor for use in minimally invasive surgery. It provides improved access and visualisation of surgical sites.

The device has a slim line shaft that can be inserted through a 5mm port. Then, once inside the abdomen, the PretzelFlex™ forms into the innovative pretzel shape.

Key Strengths of PretzelFlex™

- + Impressive strength and stability to support large organs and tissue, ideal for obesity surgery
- + Easily inserted into YelloPort+plus™ 5mm laparoscopic port access system
- + Designed and precision-machined 'in house' to maintain rigidity and shape
- + Smooth slim-line shaft designed to ensure tissue is not penetrated or damaged
- + Innovative segment technology ensures minimal space is required inside the abdomen to form into pretzel shape
- + Ideally suited to be used with FastClamp™ endoscopic clamping system



For further detail about PretzelFlex™ please contact:

sales@surginno.co.uk

+44(0)113 230 7597